

Anti-TJP1 antibody (18-299 aa) (STJ11106811)
STJ11106811

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	7082
Gene Symbol	TJP1
Uniprot ID	Z01_HUMAN
Immunogen	
Immunogen Region	18-299 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 18-299 of human TJP1 (NP_003248.3).
Immunogen Sequence	IWEQHTVTTLHRAPGFGFGIA ISGGRDNPHFQSGETSIVIS DVLKGGPAEGQLQENDRVAM VNGVSMDNVEHAFVQQLRK SGKNAKITIRRKKKVIPVS RPDPEPVSDNEEDSYDEEIH DPRSGRSGVNNRSEKIWPR DRASRERSLSPRSDRRSVA SSQPAKPTKVTLVKSRKNEE YGLRLASHIFVKEISQDSL AARDGNIQEGDVVLKINGTVT ENMSLTDAKTLIERSKGKL