

Anti-TMPO antibody (1-410 aa) (STJ11106812)

STJ11106812

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	7112
Gene Symbol	TMPO
Uniprot ID	LAP2B_HUMAN LAP2A_HUMAN
Immunogen	
Immunogen Region	1-410 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-410 of human TMPO (NP_001027454.1).
Immunogen Sequence	MPEFLEDPSVLTKDKLKSEL VANNVTLPAGEQRKDVYVQL YLQHLTARNRPPLPAGTNSK GPPDFSSDEEREPTPVLGSG AAAAGRSRAAVGRKATKKTD KPRQEDKDDLDTLTNEDL LDQLVKYGVNPGPIVGTTRK LYEKLLKLREQGTESRSST PLTISSSAENTRQNGSNDSDRYSDNEEDSKIELKLEKRE PLKGRAKTPVTLKQRRVEHN QSYSQAGITETEWTSGSSK