

Anti-TTC3 antibody (1400-1700 aa) (STJ11106813)

STJ11106813

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	7267
Gene Symbol	TTC3
Uniprot ID	TTC3_HUMAN
Immunogen	
Immunogen Region	1400-1700 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1400-1700 of human TTC3 (NP_003307.3).
Immunogen Sequence	IASETQILEGSLGISVKSHC STGDAHTVLSESNRNDEHCG NSNNKCEVIPESTSAVTNIP HVQMVAIQVSWNIIHQEVNT EPYNPFEEERQGEISRIKEH QVLQDQLQEVYENYEQIKLK GLEETRDLEEKLRHLEENK ISKTELDWFLQDLEREIKKW QQEKKEIQERLKLKKKKIKK VSNASEMYTQKNDGKEKEHE LHLQSLSEISNTLTNEKMKI EEYIKKGKEDYEESHQRAV