

Anti-PEA15 antibody (1-130 aa) (STJ11106818)

STJ11106818

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	8682
Gene Symbol	PEA15
Uniprot ID	PEA15_HUMAN
Immunogen	
Immunogen Region	1-130 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-130 of human PEA15 (NP_003759.1).
Immunogen Sequence	MAEYGTLLQDLTNNITLEDL EQLKSACKEDIPSEKSEEIT TGSAWFSFLESHNKLDKDNL SYIEHIFEISRRPDLLTMVV DYRTRVLKISEEDELDTKLT RIPSARKYKDIIQPSEEEI IKLAPPPKKA