

Anti-IFNAR2 antibody (370-475 aa) (STJ11106823)

STJ11106823

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:100-1:500 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	3455
Gene Symbol	IFNAR2
Uniprot ID	INAR2_HUMAN
Immunogen	
Region	370-475 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 370-475 of human IFNAR2 (NP_997468.1).
Immunogen Sequence	LPEVDVELPTMPKDSPQQLLE LLSGPCERRKSPLQDPFEEE DYSSTEGSGGRITFNVDLNS VFLRVLDDEDSDDLEAPLML SSHLEEMVDPEDPDNVQSNH LLASGE



Western blot analysis of various lysates using IFNAR2 Rabbit pAb (STJ11106823) at 1:500 dilution.
Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (STJUS000856) at 1:10000 dilution.
Lysates/proteins: 25 Mu g per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit
Exposure time: 30s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081