

Anti-TLK2 antibody (1-220 aa) (STJ11106831)

STJ11106831

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	11011
Gene Symbol	TLK2
Uniprot ID	TLK2_HUMAN
Immunogen	
Immunogen Region	1-220 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human TLK2 (NP_006843.2). MMEELHSLDPRRQELLEARF TGVGVSKGPLNSESNSQLC SVGSLSKKEVETPEKKQNDQ RNRKRKAEPYETSQKGTPR GKHISDYFEFAGGSAPGTSP GRSVPPVARSSPQHSLSNPL PRRVEQPLYGLDGSAAKEAT EEQSALPTLMSVMLAKPRLD TEQLAQRGAGLCFTFVSAQQ NSPSTGSGNTEHSCSSQKQ ISIQHRQTQSDLTIEKISAL
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.

St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081