

Anti-SSX1 antibody (1-188 aa) (STJ11106839)

STJ11106839

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB: 1:500-1:2000
Range	ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	
Storage	
Instruction	

TARGET INFORMATION

Gene ID	6756
Gene Symbol	SSX1
Uniprot ID	SSX1_HUMAN
Immunogen	
Immunogen Region	1-188 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-188 of human SSX1 (NP_001265620.1).
Immunogen Sequence	MNGDDTFAKRPRDDAKASEK RSKAFDDIATYFSKKEWKKM KYSEKISYVYMKRNYKAMTK LGFKVTLPPFMCNKQATDFQ GNDFDNDHNRRIQVEHPQMT FGRLHRIIPKIMPKPAEDE NDSKGVSEASGPQNDGKQLH PPGKANISEKINKRSGPKRG KHAWTHRLRERKQLVIYEEI SDPEEDDE