

Anti-Sorbitol Dehydrogenase antibody (100-200 aa) (STJ11106843)

STJ11106843

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	6652
Gene Symbol	SORD
Uniprot ID	DHSO_HUMAN
Immunogen	
Immunogen Region	100-200 aa
Specificity	A synthetic peptide corresponding to a sequence within amino acids 100-200 of human Sorbitol Dehydrogenase (NP_003095.2).
Immunogen Sequence	RENDEFCKMGRYNLSPSIF CATPPDDGNLCRFYKHNAAF CYKLPDNVTFEEGALIEPLS VGIHACRRGGVTLGHKVLVC GAGPIGMVTLLVAKAMGAAQ V