

Anti-IFIH1 antibody (1-205 aa) (STJ11106845)

STJ11106845

GENERAL INFORMATION



Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB: 1:500-1:2000
Range	IHC-P: 1:50-1:200
	ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	64135
Gene Symbol	IFIH1
Uniprot ID	IFIH1_HUMAN
Immunogen	
Immunogen Region	1-205 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-205 of human IFIH1 (NP_071451.2).
Immunogen Sequence	MSNGYSTDENFRYLISCFRA RVKMYIQVEPVLDTFLPA EVKEQIQRVTATSGNMQAVE LLLSTLEKGVWHLGWTFRELV EALRRTGSPLAARYMNPETL DLPSPSFENAHDEYLQLNL LQPTLVDKLLVRDVLDCME EELLTIEDRNRIAAEENNGN ESGVRELLKRIVQKENWFSA FLNVLRTQGNLVEQLTGS DCSES



Immunohistochemistry analysis of paraffin-embedded Human colon carcinoma using IFIH1 Rabbit pAb (STJ11106845) at dilution of 1:200 (40x lens). Microwave antigen retrieval performed with 0.01M PBS Buffer (pH 7.2) prior to immunohistochemistry staining.



Immunohistochemistry analysis of paraffin-embedded Human esophageal cancer using IFIH1 Rabbit pAb (STJ11106845) at dilution of 1:200 (40x lens). Microwave antigen retrieval performed with 0.01M PBS Buffer (pH 7.2) prior to immunohistochemistry staining.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.

St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081