

Anti-MEX3A antibody (261-340 aa) (STJ11106846)

STJ11106846

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	1576
Gene Symbol	CYP3A4
Uniprot ID	CP3A4_HUMAN
Immunogen	
Immunogen Region	261-340 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 261-340 of human MEX3A (NP_001087194.1).
Immunogen Sequence	SRDRDPVFEITGAPGNVERA REEIETHIAVRTGKILEYNN ENDFLAGSPDAAIDSRYSDA WRVHQPGCKPLSTFRQNSLG