

Anti-MBD1 antibody (386-655 aa) (STJ11106852)

STJ11106852

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	4152
Gene Symbol	MBD1
Uniprot ID	MBD1_HUMAN
Immunogen	
Immunogen Region	386-655 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 386-655 of human MBD1 (NP_001191065.1).
Immunogen Sequence	LPSVWSESEDGAGSPPPYRR RKRPSASARRHHLGPTLKPTL ATRTAQPDHTQAPTKQEAGG GFVLPPPGLDLVFLREGASS PVQVPGPVAASTEALLQEAQ CSGLSWVVVALPVKQEKADT QDEWTPGTAVLTSPVLVPGC PSKAVDPGLPSVKQEPPDPE EDKEENKDDSAKLAPEEEA GGAGTPVITEIFSLGGTRFR DTA VWLPSLQGRHSGREDGC KVVETEDTVEPTSTSWNPR