

Anti-CDKN1A/p21 antibody [SRM] (STJ11106855)

STJ11106855

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:5000-1:20000 IF/ICC:1:1000-1:4000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	1026 12575
Gene Symbol	CDKN1A Cdkn1a
Uniprot ID	CDN1A_HUMAN CDN1A_MOUSE
Immunogen	
Immunogen Region	
Specificity	Recombinant protein of human/mouse CDKN1A/p21.
Immunogen Sequence	