

Anti-Fatty Acid Synthase antibody (2212-2511 aa) [SRM] (STJ11106858)

STJ11106858

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:2000-1:6000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	2194
Gene Symbol	FASN
Uniprot ID	FAS_HUMAN
Immunogen	
Immunogen Region	2212-2511 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 2212-2511 of human Fatty Acid Synthase (FASN) Rabbit PolymAb® (NP_004095.4).
Immunogen Sequence	QQQTQLNLRSLLVNPEGPTL MRLNSVQSSERPLFLVHPLE GSTTVFHSLSRLSIPTYGL QCTRAAPLDSIHSLSLAAYYID CIRQVQPEGPYRVAGYSYGA CVAFEMCSQLQAQQSPAPTH NSLFLFDGSPTYVLAYTQSY RAKLTPGCEAEAEATEAICFF VQQFTDMEHNRVLEALLPLK GLEERVAAAVDLIKSHQGL DRQELSFAARSFYKLRRAE QYTPKAKYHGNVMLLRAKT