

## Anti-AIF1/IBA1 antibody (1-147 aa) [SRM] (STJ11106865)

STJ11106865

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	
<b>Applications</b>	WB/IHC-P/ELISA
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human/Mouse/Rat

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	SRM
<b>Concentration</b>	Lot specific
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Affinity purification
<b>Dilution</b>	WB:1:500-1:1000
<b>Range</b>	IHC-P:1:500-1:1000
	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
<b>Formulation</b>	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
<b>Isotype</b>	IgG
<b>Storage</b>	
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	199
<b>Gene Symbol</b>	AIF1
<b>Uniprot ID</b>	AIF1_HUMAN
<b>Immunogen</b>	
<b>Immunogen Region</b>	1-147 aa
<b>Specificity</b>	Recombinant fusion protein containing a sequence corresponding to amino acids 1-147 of human AIF1/IBA1 Rabbit PolymAb® (NP_001614.3).
<b>Immunogen Sequence</b>	MSQTRDLQGGKAFGLLKAQQ EERLDEINKQFLDDPKYSSD EDLPSKLEGFKEKYMEFDLN GNGDIDIMSLKRMLEKLGVP KTHLELKKLIGEVSSGSGET FSYPDFLRMMLGKRSAILKM ILMYEEKAREKEKPTGPPAK KAISELP