

Anti-AIF1/IBA1 antibody (1-147 aa) [SRM] (STJ11106866)

STJ11106866

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	IHC-P:1:1000-1:5000 IF/ICC:1:50-1:200 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	199
Gene Symbol	AIF1
Uniprot ID	AIF1_HUMAN
Immunogen	
Immunogen Region	1-147 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-147 of human AIF1/IBA1 Rabbit PolymAb® (NP_001614.3).
Immunogen Sequence	MSQTRDLQGGKAFGLLKAQQ EERLDEINKQFLDDPKYSSD EDLPSKLEGFKEKYMFDLN GNGDIDIMSLKRMLEKLGVP KTHLELKKLIGEVSSGSGET FSYPDFLRMMMLGKRSAILKM ILMYEEKAREKEKPTGPPAK KAISELP