

Anti-ICAM-1/CD54 antibody [SRM] (STJ11106871)

STJ11106871

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:2000-1:10000 IHC-P:1:100-1:500 IF/ICC:1:100-1:500 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	15894 3383
Gene Symbol	Icam1 ICAM1
Uniprot ID	ICAM1_MOUSE ICAM1_HUMAN
Immunogen	
Immunogen Region	
Specificity	Recombinant protein of human/mouse ICAM-1/CD54.
Immunogen Sequence	