

Anti-Xanthine Oxidase antibody (1000-1100 aa) [SRM] (STJ11106879)

STJ11106879

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 IHC-P:1:50-1:200 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	7498
Gene Symbol	XDH
Uniprot ID	XDH_HUMAN
Immunogen	
Immunogen Region	1000-1100 aa
Specificity	A synthetic peptide corresponding to a sequence within amino acids 1000-1100 of human Xanthine Oxidase (XDH). (NP_000370.2).
Immunogen Sequence	CIIPTKFGISFTVPFLNQAG ALLHVYTDGSVLLTHGGTEM GQGLHTKMVQVASRALKIPT SKIYISETSTNTVPNTSPTA ASVSADLNGQAVYACQTIL K