

Anti-FOXC2 antibody (402-501 aa) (STJ11106909)

STJ11106909

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	2303
Gene Symbol	FOXC2
Uniprot ID	FOXC2_HUMAN
Immunogen	402-501 aa
Immunogen Region	402-501 aa
Specificity	A synthetic peptide corresponding to a sequence within amino acids 402-501 of human FOXC2 (NP_005242.1).
Immunogen Sequence	