

Anti-LMX1A antibody (1-93 aa) (STJ11106911)
STJ11106911**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	4009
Gene Symbol	LMX1A
Uniprot ID	LMX1A_HUMAN
Immunogen	1-93 aa
Immunogen Region	1-93 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-93 of human LMX1A (NP_796372.1).
Immunogen Sequence	