

Anti-ZmSNF2 antibody (901-1089 aa) (STJ11106916)
STJ11106916

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Maize

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen	901-1089 aa
Region	
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 901-1089 of maize ZmSNF2 (XP_020405944.1).
Immunogen	FSQFLEHHVIEQQLTIAGI IYSGMYSPMLASKVGLSES EKRFWEDLDGMVRAVPTNE KLSVGDLNGHVGSTNAGYEL
Sequence	AHGGFGTHRDQAKIARMKW WKLKGETFEVFRERVFVEGI WFEEDANMWVKMATCIRKVV SKVFGVTKGSSGESKDTWWW TEDVQKAIKEKKCYRTLFI DNSTVNIER



Western blot analysis of lysates from wild type (WT) and Maize transfected with ZmSNF2 using ZmSNF2 Rabbit pAb (STJ11106916) at 1:900 dilution. A Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (STJ5000856) at 1:10000 dilution. Lysates/proteins: 25A Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL enhanced kit. Exposure time: A 180s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081