

Anti-Rap1A antibody (1-184 aa) (STJ11106920)

STJ11106920

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	ELISA
Host/Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	ELISA: Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	5906
Gene Symbol	RAP1A
Uniprot ID	RAP1A_HUMAN
Immunogen	
Immunogen Region	1-184 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-184 of human Rap1A (NP_002875.1).
Immunogen Sequence	MREYKLVVLGSGGVGKSALT VQFVQGIFVEKYDPTIEDSY RKQVEVDCQQCMLEILDTAG TEQFTAMRDLYMKNQGQGFAL VYSITAQSTFNDLQDLREQI LRVKDTEDVPMILVGNKCDL EDERVVGKEQGQNLARQWCN CAFLESSAKSKINVNEIFYD LVRQINRKTPVEKKPKKKKS CLLL