

Anti-IL7 antibody (26-177 aa) [SRM] (STJ11106939)

STJ11106939

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	IHC-P:1:50-1:200
	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	3574
Gene Symbol	IL7
Uniprot ID	IL7_HUMAN
Immunogen	
Immunogen Region	26-177 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 26-177 of human IL7. (NP_000871.1).
Immunogen Sequence	DCDIEGKDGKQYESVLMVSI DQLLDSMKEIGSNCLNNEFN FFKRHICDANKEGMFLFRAA RKL RQFLKMNSTGDFDLHLL KVSEGTITILLNCTGQVKGRK PAALGEAQPTKSLEENKSLK EQKKLNDLCFLKRLQEIKT CWNKILMGTKEH