

Anti-CDK4 antibody (204-303 aa) [SRM] (STJ11106945)

STJ11106945

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/IF/ICC/IP/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:1000-1:12000 IHC-P:1:200-1:800 IF/ICC:1:200-1:800 IP:0.5 Mu g-4 Mu g antibody for 200 Mu g-400 Mu g extracts of whole cells ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your sp
Formulation	PBS with 0.09% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	1019
Gene Symbol	CDK4
Uniprot ID	CDK4_HUMAN
Immunogen	
Immunogen	204-303 aa
Region	
Specificity	A synthetic peptide corresponding to a sequence within amino acids 204-303 of human CDK4 (NP_000066.1).
Immunogen	FAEMFRRKPLFCGNSEADQL GKIFDLIGLPPEDDWPRDVS LPRGAFPPRGPRPVQSVVPE MEESGAQLLLEMLTFNPHKR
Sequence	ISAFRALQHSYLHKDEGNPE