

Anti-MET antibody (960-1161 of aa) [SRM] (STJ11106977)

STJ11106977

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB: 1:500-1:1000
Range	IF/CC: 1:50-1:200
	ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	4233
Gene Symbol	MET
Uniprot ID	MET_HUMAN
Immunogen	
Immunogen Region	960-1161 of aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 960-1161 of human MET (NP_000236.2).
Immunogen Sequence	QIKDLGSELVRYDARVHTPH LDRLVARSVSPTTEMVSNE SVDYRATFPEDQFPNSSQNG SCRQVQYPLTDMSPILTSGD SDISSPLLQNTVHIDLSALN PELVQAVQHVVIGPSSLIVH FNEVIGRGHFGCVYHGTLTD NDGKKIHCAVKSLNRITDIG EVSQFLTEGIIMKDFSHPNV LSLLGICLRSEGSPVLVLPY MKHGDLRNFIRNETHNPTVK DLIGFGLQVAKGMKYLASK



Western blot analysis of lysates from Mouse kidney using MET Rabbit mAb (STJ11106977) at 1:1000 dilution.
 Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (STJ5000856) at 1:10000 dilution.
 Lysates/proteins: 25 µg per lane.
 Blocking buffer: 3% nonfat dry milk in TBST.
 Detection: ECL Basic Kit
 Exposure time: 10s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081