

Anti-Smad9 antibody (209-226 aa) (STJ11106980)
STJ11106980**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	
Applications	WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	4093
Gene Symbol	SMAD9
Uniprot ID	SMAD9_HUMAN
Immunogen	209-226 aa
Immunogen Region	209-226 aa
Specificity	A synthetic peptide corresponding to a sequence within amino acids 209-226 of human Smad9 (NP_005896.1).
Immunogen Sequence	