

Anti-HADHA antibody (545-763 aa) [SRM] (STJ11106984)

STJ11106984

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/IF/ICC/IP/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	IHC-P:1:50-1:200 IF/ICC:1:50-1:200 IP:0.5 Mu g-4 Mu g antibody for 200 Mu g-400 Mu g extracts of whole cells ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specif
Formulation	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	3030
Gene Symbol	HADHA
Uniprot ID	ECHA_HUMAN
Immunogen	
Immunogen Region	545-763 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 545-763 of human HADHA (NP_000173.2).
Immunogen Sequence	FYTRCLAPMMSEVIRILQE GVDPKKLDLSLTSFGFPVGA ATLVDEVGVDVAKHVAEDLG KVFGERFGGGNPELLTQMVS KGFLGRKSGKGFYIYQEGVK RKDLNSDMDSILASLKLPPK SEVSSDEDIQFRLVTRFVNE AVMCLQEGILATPAEGDIGA VFGLGFPPCLGGPFRFVDLY GAQKIVDRLKKYEAAYGKQF TPCQLLADHANSPNKKFYQ