

Anti-Acetyl-SDHA-K335 antibody (STJ11106990)

STJ11106990

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	
Host/Source	Rabbit
Reactivity	

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	6389
Gene Symbol	SDHA
Uniprot ID	SDHA_HUMAN
Immunogen	
Immunogen Region	
Specificity	Recombinant protein of human SDHA.
Immunogen Sequence	