

Anti-NCX1/SLC8A1 antibody (STJ11107000)

STJ11107000

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	6546
Gene Symbol	SLC8A1
Uniprot ID	NAC1_HUMAN
Immunogen	
Immunogen Region	
Specificity	A synthesized peptide derived from human NCX1/SLC8A1.
Immunogen Sequence	