

Anti-RPS6 antibody (1-249 aa) [SMM] (STJ11107006)
STJ11107006

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SMM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	6194
Gene Symbol	RPS6
Uniprot ID	RS6_HUMAN
Immunogen	
Immunogen Region	1-249 aa
Specificity	A synthetic peptide corresponding to a sequence within amino acids 1-249 of human S6 Ribosomal Protein (NP_001001.2). MKLNISFPATGCQKLIIEVDD ERKLRTFYEKRMATEVAADA LGEEWKGYVVRISGGNDKQG FPMKQGVLTGHRVRLLSKG
Immunogen Sequence	HSCYRPRRTGERKRKSVRGV IDANLSVLNLVIVKKGEKD IPGLTDTTPRRLGPKRASR IRKLFNLKSKEDDVRQYVVRK PLNKEGKKPRTKAPKIQRLV TPRVLQHKRRRIALKKQRTK KNKEEAAYAKLLAKRMKEA KEKRQEIQAKRRRLSSLRA



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081