

Anti-LAMP1/CD107a antibody [SRM] (STJ11107019)

STJ11107019

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:25000-1:100000 IHC-P:1:50-1:200 IF/ICC:1:50-1:200 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	16783 3916
Gene Symbol	Lamp1 LAMP1
Uniprot ID	LAMP1_MOUSE LAMP1_HUMAN
Immunogen	
Immunogen Region	
Specificity	Recombinant protein of human/mouse CD107a.
Immunogen Sequence	