

Anti-CD324/E-Cadherin antibody (157-369 aa) [SRM] (STJ11107021)

STJ11107021

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:1000-1:4000
Range	IHC-P:1:700-1:7000 IF/ICC:1:200-1:800 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	12550
Gene Symbol	Cdh1
Uniprot ID	CADH1_MOUSE
Immunogen	
Immunogen Region	157-369 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 157-369 of mouse CD324/E-Cadherin (NP_033994.1).
Immunogen Sequence	DWVIPPISCPENEKGEFPKN LVQIKSNRDKETKVFSITG QGADKPPVGVFIERETGWL KVTQPLDREAIKYILYSHA VSSNGEAVEDPMEIVITVD QNDNRPEFTQEVFEGSVAEG AVPGTSVMKVSATDADDDVN TYNAAIAYTIVSQDPELPHK NMFTVNRDTGVISVLTSGLD RESYPTYTLVQAAADLQEGE LSTTAKAVITVKD