

Anti-BIK antibody (1-130 aa) (STJ11107023)

STJ11107023

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000
	ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	638
Gene Symbol	BIK
Uniprot ID	BIK_HUMAN
Immunogen	
Immunogen Region	1-130 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-130 of human BIK (NP_001188.1).
Immunogen Sequence	MSEVRPLSRDILMETLLYEQ LLEPPTMEVLGMTDSEEDLD PMEDFDSLECMEGSDALALR LACIGDEMDVSLRAPRLAQL SEVAMHSLGLAFIYDQTEDI RDVLRSFMDGFTTLKENIMR FWRSPNPGSW



Western blot analysis of lysates from Raji cells, using BIK Rabbit pAb (STJ11107023).
Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (STJ5000656) at 1:10000 dilution.
Lysates/proteins: 25 µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081