

Anti-CD294 antibody (1-32 aa) (STJ11107024)
STJ11107024**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	
Applications	WB
Host/Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	14764
Gene Symbol	Ptgdr2
Uniprot ID	PD2R2_MOUSE
Immunogen	
Immunogen Region	1-32 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-32 of mouse CD294 (NP_034092.1).
Immunogen Sequence	MANVTLKPLCPLLEEMVQLP NHSNSSLRYIDH