

Anti-KPNA6 antibody (69-169 aa) [SRM] (STJ11107033)

STJ11107033

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:3000-1:20000 IHC-P:1:1000-1:5000 IF/ICC:1:50-1:200 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	23633
Gene Symbol	KPNA6
Uniprot ID	IMA7_HUMAN
Immunogen	
Immunogen Region	69-169 aa
Specificity	A synthetic peptide corresponding to a sequence within amino acids 69-169 of human KPNA6 (NP_036448.1).
Immunogen Sequence	DSYVSSTTGESVITREMVEM LFSDDSDLQLATTQKFRKLL SKEPSPPIDEINTPRVVD R FVEFLKRNENCTLQFEAAWA LTNIASGTSQQTKIVIEAGA V