

Anti-TGF beta 1 antibody [SRM] (STJ11107037)

STJ11107037

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:1000-1:5000 IHC-P:1:200-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	7040
Gene Symbol	TGFB1
Uniprot ID	TGFB1_HUMAN
Immunogen	
Immunogen Region	
Specificity	Recombinant protein of human TGF beta 1.
Immunogen Sequence	