

Anti-MBNL1 antibody (280-380 aa) (STJ11107046)
STJ11107046

GENERAL INFORMATION



Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	4154
Gene Symbol	MBNL1
Uniprot ID	MBNL1_HUMAN
Immunogen	
Immunogen	280-380 aa
Region	
Specificity	A synthetic peptide corresponding to a sequence within amino acids 280-380 of human MBNL1 (NP_997176.1).
Immunogen	PLEATFDLGIPQAVLPPLPK RPALEKTNNGATAVNTGIFQ YQQALANMQLQQHTAFLPPV PMVHGATPATVSAATTSATS
Sequence	VPFAATATANQIPIISAEHL T



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081