

Anti-AS160 antibody (190-450 aa) (STJ11107051)

STJ11107051

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	IP/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	IP:0.5 Mu g-4 Mu g antibody for 200 Mu g-400 Mu g extracts of whole cells
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	9882
Gene Symbol	TBC1D4
Uniprot ID	TBCD4_HUMAN
Immunogen	
Immunogen Region	190-450 aa
Specificity	Recombinant Protein corresponding to a sequence within amino acids 190-450 of human AS160 (NP_055647.2).
Immunogen Sequence	KDNEDAFYNSQKFEVLYCGK VTVTHKKAPSSLIDDCMEKF SLHEQQRLKIQGEQRGPDGP EDLADLEVVPGPSGDCLPE EADGTDTHLGLPAGASQPAL TSSRVCFPERILED SGFDEQ QEFRSRCSSVTGVQRRVHEG SQKSQPRRRHASAPSHVQPS DSEKNRTMLFQVGRFEINLI SPDTKSVVLEKNFKDISSCS QGIKHVDHFGFICRESPEPG LSQYICYVFCASESLVDE