

Anti-MMP14 antibody (21-538 aa) (STJ11107056)

STJ11107056

GENERAL INFORMATION



Product Type	Primary antibodies
Short Description	
Applications	WB/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB: 1:500-1:1000 IF/ICC: 1:50-1:200 ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	4323
Gene Symbol	MMP14
Uniprot ID	MMP14_HUMAN
Immunogen	
Immunogen Region	21-538 aa
Specificity	Recombinant Protein corresponding to a sequence within amino acids 21-538 of human MMP14 (NP_004986.1).
Immunogen Sequence	ALASLGSAQSSSFPEAWLQ QYGYLPPGDLRHTQRSPQS LSAAIAAMQKFYGLQVTGKA DADTMKAMRRPRCGVPDKFG AEIKANVRRKRYAIQGLKWQ HNEITFCIQNYTPKVG EYAT YEAIKAFRWVESATPLRFR EYPYAYIREGHEKQADIMIF FAEGFHGDSTPFDGEGGFLA HAYFPGPNIGGDTHFDSAEP WTVRNEDLNGNDIFLAVHE LGHALGLEHSSDPSAIMAP



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081