

Anti-SDPR/Cavin-2 antibody (326-425 aa) [SRM] (STJ11107058)

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:5000-1:30000 IHC-P:1:4000-1:40000 IF/ICC:1:3200-1:12800 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	8436
Gene Symbol	CAVIN2
Uniprot ID	CAVN2_HUMAN
Immunogen	
Immunogen	326-425 aa
Region	
Specificity	A synthetic peptide corresponding to a sequence within amino acids 326-425 of human SDPR/Cavin-2 (NP_004648.1).
Immunogen	EDLPSSSEQMPNDQEEESFAE GHSEASLASALVEGEIAEEA AEKATSRGSNSGMSNIDLT IVEDEEEESVALEQAQKVRY
Sequence	EGSYALTSEEAEERSDGDVPQ PAVLQVHQT