

Anti-Cyclin D2 antibody (181-280 aa) (STJ11107063)

STJ11107063

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IF/ICC:1:50-1:200 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	894
Gene Symbol	CCND2
Uniprot ID	CCND2_HUMAN
Immunogen	
Immunogen Region	181-280 aa
Specificity	A synthetic peptide corresponding to a sequence within amino acids 181-280 of human Cyclin D2 (NP_001750.1).
Immunogen Sequence	AQTFIALCATDFKFAMYPPS MIATGSVGAACGLQQDEEV SSLTCDALTELLAKITNTDV DCLKACQEQIEAVLLNSLQQ YRQDQRDGSKSELDQAST



Immunofluorescence analysis of paraffin-embedded mouse testis tissue using Anti-Cyclin D2 Rabbit pAb (STJ11107063) at a dilution of 1:100 (40x lens). Secondary antibody: Cy3-conjugated Goat anti-Rabbit IgG (H+L) (STJS001166) at 1:500 dilution. A Blue: DAPI for nuclear staining. Perform high pressure antigen retrieval with 0.01 M citrate buffer (pH 6.0) prior to IF staining.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081