

Anti-TIGIT antibody (22-141 aa) (STJ11107075)

STJ11107075

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000
	ELISA: Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	201633
Gene Symbol	TIGIT
Uniprot ID	TIGIT_HUMAN
Immunogen	
Immunogen Region	22-141 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 22-141 of human TIGIT (NP_776160.2).
Immunogen Sequence	MTGTIETGNISAEKGGSI LQCHLSSTTAQVTQVNWEQQ DQLLAICNADLGWHISPSFK DRVAPGPGGLTLQSLTVND TGEYFCIYHTYPDGTYTGR I FLEVLESSVAEHGARFQIP



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081