

Anti-ABCC5/MRP5 antibody (1101-1200 aa) (STJ11107079)

STJ11107079

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IF/ICC:1:50-1:200 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	10057
Gene Symbol	ABCC5
Uniprot ID	MRP5_HUMAN
Immunogen	
Immunogen Region	1101-1200 aa
Specificity	A synthetic peptide corresponding to a sequence within amino acids 1101-1200 of human ABCC5/MRP5 (NP_005679.2).
Immunogen Sequence	RDLISIALITTTGLMIVLM HGQIPPAYAGLAISYAVQLT GLFQFTVRLASETEARFTSV ERINHYIKTLSLEAPARIKN KAPSPDWPQEGEVTFENAEM



Immunofluorescence analysis of A paraffin-embedded Human kidney tissue using A ABCC5/MRP5 Rabbit pAb (STJ11107079) at A dilution of A 1:200 (40x lens). Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) (STJ5001166) at 1:500 dilution. A Blue: DAPI for nuclear staining. Perform high pressure antigen retrieval with 0.01 M citrate buffer (pH 6.0) prior to IF staining.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081