

Anti-Laminin beta 1 antibody (40-270 aa) (STJ11107089)
STJ11107089

GENERAL INFORMATION



Product Type	Primary antibodies
Short Description	
Applications	WB/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:1000-1:5000
Range	IF/CC:1:100-1:500 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	3912
Gene Symbol	LAMB1
Uniprot ID	LAMB1_HUMAN
Immunogen	
Immunogen Region	40-270 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 40-270 of human Laminin beta 1 (NP_002282.2).
Immunogen Sequence	GDLLIGRAQKLSVTSTCGLH KPEPYCIVSHLQEDKKCFIC NSQDPYHETLNPDSHLIENV VTTFAPNRLKIWWQSENGVE NVTIQLDLEAEFHFTHLIMT FKTFRPAAMLIERSSDFGKT WGVYRYFAYDCEASFPGIST GPMKKVDDIICDSRYSDIEP STEGEVIFRALDPAFKIEDP YSPRIQNLLKITNLRKFVK LHTLGDNLLDSRMEIREKYY YAVYDMVVRGN



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081