

Anti-Marchf2 pAb antibody (147-246 aa) (STJ11107093)

STJ11107093

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	224703
Gene Symbol	Marchf2
Uniprot ID	MARH2_MOUSE
Immunogen	
Immunogen Region	147-246 aa
Specificity	A synthetic peptide corresponding to a sequence within amino acids 147-246 of mouse Marchf2 (NP_001239409.1).
Immunogen Sequence	FITPLAAISGWLCLRG AQDH LRLHSRLEAVGLIALTIALF TIYVLWTLVSFRYHCQLYSE WRKTNQKVRLKIREADGSED PHHSLLATGLLKKVAEETPV