

Anti-CD3/CD28-bispecific antibody [SM] (STJ11107096)
STJ11107096**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	
Applications	FC
Host/Source	
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	FC:0.2 Mu g per 10 ⁶ cells in 100 Mu l volume
Formulation	PBS, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	940 916
Gene Symbol	CD28 CD3E
Uniprot ID	CD28_HUMAN CD3E_HUMAN
Immunogen	
Immunogen Region	
Specificity	Recombinant protein of human CD3/CD28.
Immunogen Sequence	