

Anti-CD62L/L-Selectin antibody (39-332 aa) (STJ11107101)

STJ11107101

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	IHC-P/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC-P:1:500-1:1000 IF/ICC:1:100-1:500 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	6402
Gene Symbol	SELL
Uniprot ID	LYAM1_HUMAN
Immunogen	
Immunogen Region	39-332 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 39-332 of human CD62L/L-Selectin (NP_000646.3).
Immunogen Sequence	WTYHYSEKPMNWQRARRFCR DNYTDLVAIQNKAEIYLEK TLPFSRSYYWIGIRKIGGIW TWVGTNKSLEEAENWGDGE PNNKKNKEDCVEIYKRNKD AGKWND DACHKLKAALCYTA SCQPWSCSGHGECVEIINNY TCNCDVGYYGPQCQFVIQCE PLEAPELGTMDC THPLGNFS FSSQCAFSCSEGTNLTGIEE TTCGPGFNWSSPEPTCQVIQ CEPLSAPDLGIMNCSHPLA



Immunohistochemistry analysis of paraffin-embedded Human tonsil tissue using CD62L/L-Selectin Rabbit pAb (STJ11107101) at a dilution of 1:1000 (40x lens). High pressure antigen retrieval performed with 0.01M Citrate Buffer (pH 6.0) prior to immunohistochemistry staining.



Immunofluorescence analysis of Jurkat cells using CD62L/L-Selectin Rabbit pAb (STJ11107101) at a dilution of 1:400 (40x lens). Secondary antibody: Cy3-conjugated Goat anti-Rabbit IgG (H+L) (STJ3001166) at 1:500 dilution. Blue: DAPI for nuclear staining.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081