

Anti-PPARD antibody (101-200 aa) (STJ11107104)

STJ11107104

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	5467
Gene Symbol	PPARD
Uniprot ID	PPARD_HUMAN
Immunogen	
Immunogen Region	101-200 aa
Specificity	A synthetic peptide corresponding to a sequence within amino acids 101-200 of human PPARD (NP_006229.1).
Immunogen Sequence	TIRMKLEYEKSCRSCKIQKK NRNKCQYCRFQKCLALGMSH NAIRFGRMPEAEKRKLVAAGL TANEGSQYNPQVADLKAFSK HIYNAYLKNFNMTKKKKARSI



Western blot analysis of lysates from Rat brain using PPARD Rabbit pAb (STJ11107104) at 1:1000 dilution incubated overnight at 4°C.
Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (STJ5000856) at 1:10000 dilution.
Lysates/proteins: 25 Mu g per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081