

Anti-IL12 antibody (23-328 aa) (STJ11107105)
STJ11107105

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	3593 3592
Gene Symbol	IL12B IL12A
Uniprot ID	IL12B_HUMAN IL12A_HUMAN
Immunogen	
Immunogen Region	23-328 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 23-328 of human IL12B/23-219 of human IL12A (NP_002178.2/NP_001384921.1).
Immunogen Sequence	IWELKKDVYVVELDWYPDAP GEMVVLTCDTPEEDGITWTL DQSSEVLGSGKTLTIQVKEF GDAGQYTCHKGGEVLSHSL LLHKKEDGIWSTDILKDQKE PKNKTLRLCEAKNYSGRFTC WWLTTISTDLTFSVKSSRGSDPQGVTCGAATLSAERV DNKEYEYSEVCQEDSACPAA EESLPIEVMVDAVHKLKYEN YTSSFFIRDIIKPDPPKNLQ LKPLKNSRQVEVSWEYPT



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081