

Anti-Pan-Lactic Acid-Lysine antibody [SRM] {Magnetic Beads} (STJ11107111)

STJ11107111

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	IP
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat/Other

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Magnetic Beads
Purification	Affinity purification
Dilution Range	IP:30ul antibody (bead slurry) for 200 Mu g-400 Mu g extracts of whole cells
Formulation	PBS with 0.05% Proclin300, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	A synthetic peptide corresponding to a sequence containing Pan Lactic acid-Lysine.
Immunogen Sequence	