

Anti-NNMT antibody (1-264 aa) [SRM] (STJ11107119)

STJ11107119

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/IF/ICC/IP/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SRM
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:1000-1:5000
Range	IF/CC:1:50-1:200 IP:0.5 Mu g-4 Mu g antibody for 200 Mu g-400 Mu g extracts of whole cells ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirement
Formulation	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	4837
Gene Symbol	NNMT
Uniprot ID	NNMT_HUMAN
Immunogen	
Immunogen Region	1-264 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-264 of human NNMT (NP_006160.1).
Immunogen Sequence	MESGFTSKDITYLSHFNPRDY LEKYKFGSRHSAESQILKH LLKNLFKIFCLDGVKGDLLI DIGSGPTIYQLLSACESFKE IVTDYSDQNLQELEKWLLK EPEAFDWSPWVTYVCDLEGN RVKGPEKEEKLQAVKQVLK CDVTQSQPLGAVPLPPADCV LSTLCLDAACPDLPITYCRAL RNLGSLLKPGGFLVIMDALK SSYYMIGEQKFSSLPLGREAE VEAHVKEAGYTIEWFEVIS